

## High Resolution Chromosomes

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

### Reagents

**Acetic acid, glacial**

**Bromodeoxyuridine (BrdU)**

Sigma, Cat. B9285

**Colcemid, KaryoMAX Colcemid Solution (10 µg/ml)**

Gibco BRL, Cat.15210-016

**Fetal Bovine Serum Qualified, heat inactivated**

Gibco BRL, Cat. 16140-022, 500 ml

**L-Glutamine-200 mM, 100x**

Gibco BRL, Cat. 25030-016 (100 ml)

**Penicillin/Streptomycin 5,000 U/ml/5,000 µg/mL**

Gibco BRL, Cat. 15070-014

**Phytohaemagglutinin (PHA), HA 15**

Murex Diagnostics Ltd., Dartford, England DA1 5LR

**Methotrexate (MTX)**

Sigma, Cat. M8407

**Methanol**

**Potassium chloride**

Mallinckrodt, Cat. 6858

**RPMI Medium 1640**

Gibco BRL, Cat. 21870-050

### Preparation

#### RPMI 1640 Full Medium

Components	Amount
RPMI Medium 1640	385 ml
L-Glutamine-200 mM, 100x	5 ml
Penicillin/Streptomycin 5,000 U/ml/5,000 µg/ml	10 ml
Fetal Bovine Serum Qualified, heat inactivated	100 ml

**Hypotonic solution**

0.075 M potassium chloride (KCl) in distilled water

**Fixative**

Methanol/glacial acetic acid, 3:1(v/v)

**MTX stock**

$10^{-5}$  M in H<sub>2</sub>O

**BrdU stock**

1 mg/ml in distilled water

**Procedure**

1. Use T25 (with 5 ml media) or T75 flasks (with 20 ml media).
2. Initiate PHA-stimulated lymphocyte cultures. Incubate in upright position at 37°C.
3. At 72 h add from the MTX stock ( $10^{-5}$  M) to a final concentration of  $10^{-7}$  M; mix well and incubate an additional 17 h.
4. After 17 h centrifuge the content of the flasks, remove the supernatant, and wash the pellet twice with plain media.
5. After the second wash resuspend the pellet in RPMI 1640 20% BSA and transfer to a flask.
6. Add from the BrdU stock (1 mg/ml) to a final concentration of 25 µg/ml.
7. Incubate for 5 h 30 min at 37°C.
8. During the last 10 min. of incubation add Colcemid stock (10 µg/ml) to a final concentration of 0.06 µg/ml.
9. Centrifuge cultures for 10 min.
10. Remove supernatant and add hypotonic solution, 0.075M KCl; incubate for 10 min at 37°C.
11. Add a few drops of fresh fixative and spin for 10 min.
12. Aspirate supernatant and add fixative.

13. Resuspend pellet gently in fixative and centrifuge for 5 min. Repeat this step two more times.

14. Store pellet under fixative at -20°C until ready to prepare slides.