

## **Metaphase Preparation from Mouse/ Rat Spleen**

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

### **Reagents**

**Acetic acid, glacial**

Mallinckrodt, Cat. V193

**Antibiotic-Antimycotic, 100X**

10,000 U/ml Penicillin G sodium; 10,000 µg/ml streptomycin sulfate;  
25 µg amphotericin B Gibco BRL, Cat 15240-013

**or**

**Penicillin/ Streptomycin**

5,00 U/ml/5,000 µg/ml

Gibco BRL, Cat.15070-014

**β-Mercaptoethanol**

Sigma, Cat. 6250

**Concanavalin A, 5 mg**

Sigma-Aldrich, Cat. C-5275

**Colcemid, KaryoMAX Colcemid Solution, 10 µg/ml**

Gibco BRL, Cat. 15210-016

**Fetal Bovine Serum Qualified, heat inactivated**

Gibco BRL, Cat. 16140-022

**Homogenizer**

Thomas Scientific, Cat. 3431D7

**L-Glutamine-200 mM, 100X**

Gibco BRL, Cat. 25030-016

**Lipopolysaccharides (LPS), 25 mg**

Sigma-Aldrich, Cat. L-2637

**Methanol alcohol, Anhydrous**

Mallinckrodt, Cat. 3016

**Potassium Chloride (KCl)**

Mallinckrodt, Cat. 6858

**RPMI Medium 1640**

Gibco BRL, Cat. 21870-050

## Preparation

### Concanavalin A (5mg)

Stock solution: 5  $\mu\text{g}/\mu\text{l}$  in tissue culture grade water

Store aliquots at  $-20^{\circ}\text{C}$

### Lipopolysaccharides (LPS) (25mg)

Stock solution: dissolve 25 mg in 1 ml sterile water

Use 1:1000 dilution for a final concentration of 25  $\mu\text{g}/\text{ml}$  of culture

### 0.5% Beta- Mercaptoethanol

2 $\mu\text{l}$  100% B-M, 400  $\mu\text{l}$  sterile water

### Hypotonic Solution (0.075M)

KCl 5.6 g

Distilled water 100 ml

Pre-warm to  $37^{\circ}\text{C}$

### Fixative

Prepare fresh methanol:acetic acid, 3:1

### Medium

RPMI 1640 500 ml

Antibiotic –Antimycotic 100x 5 ml

**or**

Penicillin/Streptomycin 5 ml

L- glutamine 5 ml

Do not add FBS!!

Filter sterilize solution

## Procedure

1. Transport spleen in medium w/o FBS.
2. Place a single spleen into a homogenizer with 3 ml of prepared medium w/o FBS. Grind well.
3. Transfer 0.5 ml to one T75 Nunc flask.
4. Add 20 ml of prepared media.

5. Add the following to this flask:

|                               |        |
|-------------------------------|--------|
| FBS                           | 2.5 ml |
| Concanavalin A, 5 µg/ul       | 30 µl  |
| Lipopolysaccharides (diluted) | 25 µl  |
| β-Mercaptoethanol, 0.5%       | 30 µl  |

6. Incubate flask at 37°C for 48 hours.

7. At the end of the 48-hour incubation, add 0.25 ml Colcemid to each flask. Continue to incubate for an additional 30-60 min until cells start to divide.

8. Transfer the suspension to 50 ml tubes.

9. Centrifuge at 1,000 rpm for 10 min.

10. Remove supernatant leaving 500 µl of solution.

11. Resuspend pellet.

12. Gently add, drop-wise, 5 ml 0.075 M KCl, prewarmed to 37°C.

13. Incubate at 37°C for 20 minutes.

14. Add 1 ml of freshly prepared fixative.

15. Centrifuge at 1,000 rpm for 10 minutes.

16. Remove supernatant to 500 µl.

17. Resuspend pellet.

18. Transfer to 15 ml tubes.

19. Add fixative to 10 ml.

20. Centrifuge at 1,000 rpms. Repeat wash with fresh fixative at least 3 x.

21. Remove the supernatant one final time and drop 12 µl of the suspension onto a clean slide in the humidity of a drying chamber or Thermotron.

22. The slides are stored in a drawer at room temperature for one week, then sealed in a heat-sealable pouch (e.g. Kapak/Scotchpak) with Drierite until use.