

Preparation of Protein Extracts for 2D Gel-analysis

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Aprotinin

Sigma, Cat. A1153

Benzamidine

Sigma, Cat. B6508

Hepes, pH 8**Isopropanol (2-Propanol)****Leupeptin**

Sigma, Cat. L2884

Phenylmethylsulfonyl fluoride (PMSF)

Sigma, Cat. P7626

1X Phosphate Buffered Saline (PBS), pH 7.4

Gibco/BRL, Cat. 10010-023

Sodium fluoride (NaF)

Sigma, Cat. S1504

Sodium orthovanadate (Na₃VO₄)

Sigma, Cat. S6508

Sodium pyrophosphate (Na₄P₂O₇)

Sigma, Cat. S6422

Water, sterile

Preparation

Aprotinin [1000X]

1 mg/ml in 0.01M HEPES, pH 8

Store at -20°C

Leupeptin [1000X]

1 mg/ml in sterile water

Store at -20°C

PMSF [50X]

1.74 mg/ml in isopropanol

Store at -20°C

2-D Lysis Buffer

1X PBS	100	ml	
Na ₄ P ₂ O ₇	223	mg	f.c. [5mM]
Na ₃ VO ₄	1.8	mg	f.c. [100µM]
NaF	21	mg	f.c. [5mM]
Benzamidine	13	mg	f.c. [830µM]
Store at 4°C			

2-D Lysis Buffer + Protease Inhibitors (PIH)

Dilute fresh from Stock into desired aliquot of 2-D Lysis Buffer:

Aprotinin	1:1000
Leupeptin	1:1000
PMSF	1:50

Procedure

1. Wash cells 2 x with 1X PBS pre-chilled to 4°C.
2. Scrape cells off plate into ice cold 2-D Lysis Buffer + PIH. Volume of 2-D Lysis Buffer + PIH should be determined by the size of plate/flask in which the cells are growing; 2.5 ml is sufficient for T-150 flask.
3. Pellet at 660 x g for 3 min (about 2000 rpm in clinical centrifuge) at 4°C.
4. Resuspend pellet in 1 ml 2-D Lysis Buffer + PIH, transfer to pre-weighed eppendorf tube and pellet at 2700 x g for 5 min (about 5700 rpm in microfuge) at 4°C.
5. Aspirate off 2-D Lysis Buffer and determine weight of pellet.
6. Store at -80°C until further processing for 2-D gels.

Notes

1. If wet weight (WW) of pellet is 5-10 mg (~10 million cells) it is usually possible to obtain one acceptable 2-D gel. It is always recommended to obtain cell pellets of more than 15 mg, so it may be necessary to harvest and pool together many plates in order to obtain a sufficient pellet.

