

Amplification of Genomic DNA using Alu PCR

Reagents

DNA fragment block from gel or YAC DNA

dATP, dTTP, dGTP, dCTP

Boehringer Mannheim, Cat. 105 1440, 105 1458, 105 1466, 105 1482

Ethanol, absolute

Sigma, Cat. E7023

MgCl₂ (25 mM)

(comes with Perkin Elmer buffer)

10X PCR Buffer II (without MgCl₂)

Perkin Elmer, Cat. N808-0010

Primer CL1 (50 μM)

(5'- TCCCAAAGTGCTGGGATTACAG -3')

Primer CL2 (50 μM)

(5' -CTGCACTCCAGCCTGGG -3')

Sodium chloride (NaCl) (5 M)

Quality Biologicals, Cat. 351-036-100

Taq Polymerase: Ampli Taq (5 U/μl)

Perkin Elmer, Cat. N808-0158

Tris-HCl (1 M), pH 7.5

Quality Biologicals, Cat. 351-006-100

H₂O

Preparation

DNA

Cut desired band from agarose gel (block) or prepare DNA.

Wash gel block with 1M Tris-HCl for 3 x 45 min.

Incubate gel block with 1X PCR-buffer in order to remove EDTA.

Melt gel block in 70°C water bath for 5 min.

dNTP :Nucleotides (Total Concentration 100mM)

dATP, 100 μl for a final concentration of 25mM

dTTP, 100 μl for a final concentration of 25mM

dGTP, 100 μl for a final concentration of 25mM

dCTP, 100 μl for a final concentration of 25mM

Mix and store at -20°C.

Procedure

1. Prepare the PCR reaction mix (final volume 100 μ l)

DNA (100ng),	X μ l
10X PCR Buffer II(without MgCl ₂),	10 μ l
MgCl ₂ 25mM	6 μ l
dNTP :Nucleotides (25mM),	1 μ l
Primer CL1 (50 μ M)	0.5 μ l
Primer CL2 (50 μ M)	0.5 μ l
Sterile H ₂ O	82.5 μ l – X μ l
*Taq Polymerase	0.5 μ l

*Add Taq after the first 3 min step

2. Carry out PCR according to the following profile (Program called “DOPMS”):

<u>Step</u>	<u>Temperature (°C)</u>	<u>Minutes</u>
1	96	3
2	96	1
3	37	0.5
4	72	6
5	repeat steps 2-4, 29 times	
6	72	5
7 (cooling)	4	

Purification (steps 3-8)

3. Check 10 μ l of PCR product on a 1.2% agarose gel.
4. Precipitate the DNA using 1/25 vol. 5 M NaCl, 1/50 vol. 0.5 M MgCl₂, and 2.5 vol. absolute ethanol for at least 10 min at -80°C.
5. Centrifuge for 30 min at 13,000 rpm.
6. Pour off supernatant and wash with 400 μ l 70% ethanol.
7. Dry in the speed vac.
8. Dissolve DNA in 50 μ l sterile water.