

Detection (SKY)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Avidin-Cy5 (1.8 mg/ml)

Jackson Immuno Research Lab, Cat. 003-170-083

Bovine Serum Albumin (BSA)**DAPI**

Sigma, Cat. 18860

Formamide

Fluka BioChemika, Cat. 47671

HCl, 1 N**Mouse anti-digoxigenin (0.1 mg/ml)**

Sigma, Cat. D 8156

Sheep anti-mouse Cy5.5 (1.0mg/ml)

Amersham, Cat. RPQ 0115

SSC, 20X**Tween 20**

Sigma, Cat. P-1379

dH₂O**Antifade (1,4-phenylene-diamine)**

Preparation

50% FA/2X SSC

20X SSC 30 ml

dH₂O 120 ml

Formamide 150 ml

Adjust pH to 7.0 using 1 N HCl

Pre-warm to 45(C)**1X SSC**

20X SSC 25 ml

dH₂O 475 ml

Pre-warm to 45(C)

4X SSC/0.1%Tween 20

20X SSC 200 ml
dH₂O 799 ml
Tween 20 1 ml

Pre-warm to 45(C

Blocking Solution (3% BSA/4X SSC/0.1%Tween 20)

BSA 0.3 g
4X SSC/0.1%Tween 20 10 ml
Vortex until dissolved

Pre-warm to 37(C

Antibody Solution (1% BSA/4X SSC/0.1%Tween 20)

0.1 g BSA
10 ml 4X SSC/0.1%Tween 20

Pre-warm to 37(C

DAPI stock solution (f.c.= 0.2 mg/ml)

DAPI 2 mg
dH₂O 10 ml
Aliquot and store at -80°C

DAPI staining solution (f.c.= 80 ng/ml)

DAPI (stock solution) 40 (l
2X SSC 100 ml
Store at 4°C in a light-tight coplin jar.

Procedure

1. Carefully remove rubber cement surrounding coverslips with forceps.
2. Wash slides in 50% formamide/2X SSC for 3 x 5 min each, shaking.
3. Wash slides in 1X SSC for 3 x 5 min, shaking.
4. Dip slides in 4X SSC/0.1% Tween 20; do not let them dry.
5. Add 120 µl of Blocking Solution (3% BSA/4X SSC/0.1%Tween20) to a 24 mm x 60 mm coverslip and incubate in a moist hybridization chamber at 37(C for 30 min.
6. Dip slides in 4X SSC/0.1% Tween 20 to wash off blocking solution.
7. Spin all fluorescent dyes for 3 min at 13,000 rpm.

8. Combine the two antibodies, mouse anti-Dig and Avidin-Cy5, into the same eppendorf tube, and apply 120 μ l of antibody solution to a 24 mm x 60 mm coverslip. Each antibody should be diluted 1:200 in 1% BSA/4X SSC/0.1% Tween 20 or 4X SSC/0.1% Tween 20. Invert the slide onto the solution. Incubate the slides in a moist hybridization chamber at 37(C for 45-60 min.
9. Wash slides in 4X SSC/0.1% Tween 20, 3 x 5 min, shaking.
10. Add 120 μ l of the antibody (sheep anti-mouse Cy5.5, diluted 1:200 in 1% BSA/4X SSC/0.1% Tween 20 or 4X SSC/0.1% Tween 20). Incubate slides in a moist hybridization chamber at 37(C for 45 min.
11. Wash slides 3 x 5 min, shaking.
12. Stain slides for 5 min in the DAPI staining solution in a light-protected coplin jar at RT.
13. Wash slides with 2X SSC 3-5 min. Dehydrate slides in ethanol series of 70%, 90%, and 100% for 3 min each; air-dry slides.
14. Apply 35 μ l of antifade solution, cover each slide with a 24 mm x 60 mm coverslip, and store in a light-protected container at 4°C until slide is imaged.

Notes

1. Exposure of slides to ambient light should be minimized during all procedures.
2. Carefully remove coverslips during all procedures to minimize scratches.
3. Do not let the slide dry out between washing steps.