



Embryoid Body Formation Protocol
Adapted from:
Human Embryonic Stem Cells: Laboratory Manual
(Includes Invitrogen product information)

Michal Amit and Joseph Itskovitz-Eldor.

Department of Obstetrics and Gynecology, Rambam Medical Center and The Bruce Rappaport Faculty of Medicine, Technion-Israel Institute of Technology

The use of the procedures contained herein is the subject of Patents, Patents applications and other proprietary Intellectual Property owned by the authors and/or by third parties. This document should not be considered an express or implied license under any Intellectual Property rights.

Correspondence to Joseph Itskovitz-Eldor M.D., D.Sc., Department of Obstetrics and Gynecology, Rambam Medical Center, P.O.B. 9602, Haifa 31096, Israel.

Tel.: +972-4-854-2536; Fax: +972-4-854-2503;

E-mail: Itskovitz@rambam.health.gov.il.

Copyright © 2002 Technion – Israeli Institute of Technology. All Rights reserved. 1
The information is copyrighted in Israel and worldwide by international treaty.
Revised; October 2002.

Formation of Embryoid Bodies (EBs):

1. Remove medium from well. Add 0.5 ml splitting medium (see A below), and incubate for at least 30 minutes.
2. Add 1 ml of culture medium (see B below) and gently scrape cells with 5-ml pipette.
3. Collect cell suspension and place into conical tube.
4. Centrifuge 3 minutes at 800 rpm at a recommended temperature of 4 °C.
5. Re-suspend cells in media (see B below) using Gillson 1000 μ M tip and plate on 58 mm Petri dish.
6. Add 6 ml of medium.

Note:

If EBs attach to the dish, scrape them off gently .

A. hES cell splitting medium:

1 mg / ml [Collagenase type IV](#). Invitrogen cat. # 17104

[Dulbecco's Modified Eagle's Medium \(DMEM\)](#). Invitrogen cat #11960.

B. hES cell media:

B.1 Normal medium:

Final concentrations:

80% [Dulbecco's Modified Eagle's Medium \(DMEM\)](#). Invitrogen cat #11960
or Knockout DMEM [KO-DMEM](#). Invitrogen cat. # 10829.

20% [Fetal Bovine Serum defined \(FBSd\)](#). Invitrogen cat. # 16141

1% [Non essential amino acids](#). Invitrogen cat. # 11140

mM [L-glutamine](#) Invitrogen cat. # 21051

0.1 mM β -[Mercaptoethanol](#) Invitrogen, cat. # 21985

Preparation:

1. Pour all materials into 22 μ M filter unit, and filter.
2. Store at 4°C.

B.2 Serum free medium:

Final concentrations:

80% [KO-DMEM](#) Invitrogen cat. # 10829.

20% [Serum replacement \(SR\)](#) Invitrogen cat. # 10828.

1% [Non essential amino acids](#). Invitrogen cat. # 11140

0 mM [L-glutamine](#). Invitrogen cat. # 21051

0.1 mM β -[Mercaptoethanol](#). Invitrogen, cat. # 21985

4 ng/ml [basic Fibroblasts Growth Factor \(bFGF\)](#). Invitrogen cat. # 13256

Preparation:

1. Pour all materials into 22 μ M filter unit, and filter.
2. Store at 4°C .

May be used within two weeks of preparation.