

Hematoxylin and Eosin (H&E) staining

- Place slides containing paraffin sections in a slide holder (glass or metal)
- Deparaffinize and rehydrate sections:

3 x 3'	Xylene (<i>blot excess xylene before going into ethanol</i>)
3 x 3'	100% ethanol
1 x 3'	95% ethanol
1 x 3'	80% ethanol
1 x 5'	deionized H ₂ O
- While sections are in water, skim surface of hematoxylin with a Kimwipe to remove oxidized particles. Blot excess water from slide holder before going into hematoxylin.
- Hematoxylin staining:

1 x 3'	Hematoxylin
Rinse	deionized water
1 x 5'	<u>Tap</u> water (<i>to allow stain to develop</i>)
Dip 8-12x (fast)	Acid ethanol (<i>to destain</i>)
Rinse 2 x 1'	Tap water
Rinse 1 x 2'	Deionized water (<i>can leave overnight at this stage</i>)
- Blot excess water from slide holder before going into eosin.
- Eosin staining and dehydration:


1 x 30 sec	Eosin (<i>up to 45 seconds for an older batch of eosin</i>)
3 x 5'	95% ethanol
3 x 5'	100% ethanol (<i>blot excess ethanol before going into xylene</i>)
3 x 15'	Xylene
- You can leave slides in xylene overnight to get good clearing of any water.
- Coverslip slides using Permount (xylene based).
- Place a drop of Permount on the slide using a glass rod, taking care to leave no bubbles.
- Angle the coverslip and let fall gently onto the slide. Allow the Permount to spread beneath the coverslip, covering all the tissue.
- Dry overnight in the hood.

Reagents for H&E staining:


- **Xylene:** StatLab (Lewisville, TX) #8400
Laboratory grade, Anapath brand
- **Acid Ethanol:** 1 ml concentrated HCl + 400 ml 70% ethanol
- **Hematoxylin:** Poly Scientific (Bayshore, NY) #s212A
Harris hematoxylin with glacial acetic acid
- **Eosin:** Poly Scientific (Bayshore, NY) #s176
Eosin Phloxine stain, working
- **Permount:** Fisher Scientific #SP15-100
Histological mounting medium

Representative Images:

Mammary gland ducts (mature virgin)



20x magnification



40x magnification