

INVERSE PCR (Jill)

1. Digest 2.5 μ g of genomic DNA with restriction enzymes. Stop reaction at 65°C for 20'.
2. Fractionate on a 1% agarose gel to check DNA has been cut.
3. Set up self-ligation reaction:

40 μ l	10 x DNA Ligase buffer
1.5 μ l	NEB concentrated DNA ligase
358.5 μ l	1 μ g digested DNA plus dH ₂ O

Incubate overnight at room temperature. Do not put in more than 1 μ g of DNA as fragments may ligate to each other rather than to themselves.

4. Add 400 μ l 25:24:1 phenol:chloroform: IAA. Vortex and spin for 5' in microfuge. Remove aqueous layer to fresh tube.
5. Add 2.5 volumes ethanol and 0.1 volume 3M NaOAc pH4.5. Precipitate at -20°C overnight or at -80°C for 1 hr.
6. Spin in microfuge for 10' at full speed. Wash pellet with 70% ethanol, spin again and drain off supernate.
7. Resuspend pellet in 30 μ l dH₂O and check concentration on a gel or by OD₂₆₀
8. Set up PCR with 50ng DNA using extend taq for long products:

37.5 μ l	Ligated DNA plus dH ₂ O
5 μ l	10 X PCR buffer
4 μ l	25 mM MgCl ₂
1 μ l	10 mM dNTP's
1 μ l	Primer 1 (10 μ M)
1 μ l	Primer 2 (10 μ M)
0.5 μ l	Taq DNA Polymerase

50 μ l

TOTAL

9. Cycle:

94°C (45s)

X 1

94°C (15s), 50°C (30s), 68°C (4')

X 10

94°C (15s), 50°C (30s), 68°C (4' plus 5s per cycle)

X 20

72°C (7')

X 1

10. 25 μ l of the final PCR product should be run out on a 1% agarose gel, viewed and blotted for band verification by hybridization. The remaining 25 μ l of the PCR reaction can then be used, if positives are forthcoming, to clone and sequence the band of interest.