

Oligo Annealing

Get 64-mer at 0.05 μ M scale, desalted (no HPLC or PAGE purification necessary)

Annealing

- dissolve both oligos in 50 μ l of water
- take 2 μ l of each oligo (forward and reverse)
- add 48 μ l of annealing buffer (see below)
- incubate 4 min. at 95°C
- incubate 10 min. at 70°C
- cool down oligos on your bench (can be stored at -20°C)

Phosphorylation

- take 5 μ l of annealed oligos
- add 12 μ l of water
- add 2 μ l of T4 ligase buffer (usually containing 1mM ATP)
- add 1 μ l PNK kinase
- incubate 30 min. at 37°C
- incubate 10 min. at 70°C to heat inactivate PNK

Ligation

- take 5 μ l of phosphorylated oligos
- add 1 μ l of vector (20-100ng, CIPed)
- add 11 μ l of water
- add 2 μ l of T4 ligase buffer
- add 1 μ l of T4 ligase (1-5U)

Transform and grow bacteria as usual

Check minipreps by digestion

Sequence since some oligos may be mutated.

Annealing buffer:

- 100mM potassium acetate
- 30mM HEPES pH7.4
- 2mM magnesium acetate