OXFORD Practical Approach Series

Page 1 of 1

First strand cDNA synthesis

Michael L. Nuccio, Tzung-Fu Hsieh, and Terry L. Thomas

Department of Biology, Texas A&M University, College Station, Texas 77843-3258, U.S.A.

Equipment and reagents

- cDNA synthesis system (Gibco BRL)
- DEPC treated water
- 25 μM T_{11 or 12}VN oligonucleotide primer

Method

- 1 Mix the following reagents in a microcentrifuge tube:
 - 1 μl total RNA (2 μg)
 - 5 μ l of 5 × reverse transcriptase reaction buffer
 - 2.5 μl of 25 μM T_{11 or 12}VN primer
 - 2.5 μl of 200 μM dNTP mix
 - 11.5 μI DEPC H₂O
- 2 Mix by pipetting gently, then heat to 65 °C for 3 min, then let stand at room temperature for 3 min. Repeat this step one more time.
- 3 Allow the reaction mix to cool to room temperature then add:
 - 1.25 μI of 0.1 M DTT
 - 1.25 µl Superscript (250 U)
- 4 Mix by pipetting.
- 5 Incubate at 42 °C for 60 min.
- 6 Heat inactivate reaction by incubating at 95 °C for 5 min.

[©] Oxford University Press 23 August, 2001 All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, without the prior permission in writing of Oxford University Press, or as expressly permitted by law, or under terms agreed with the appropriate reprographics rights organization.